







## **CLINNOVA** Gynae Pro









The CLINNOVA Gynae Pro is intelligently designed for colposcopy & gynaecology procedures, offering fully electric profiling of height, backrest and two-way Trendelenburg tilt. The couch also includes a CPR facility and digital display handset, with safe working load of 285kg.



### **Standard Features**

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 285Kg (627lbs)
- Electric backrest, adjustable from 0° to +90°
- Electric height elevation, raises from 52cm to 100cm in 15 seconds
- Electric trendelenburg and reverse trendelenburg
- Choice of independent, retractable or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Includes a digital display handset as standard
- A pair of multi-adjustable leg supports and arm loops supplied as standard
- Seat section lifts for easy cleaning of framework
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

# **Optional Accessories**

9210 10 Button Foot Switch
9043 Removable leg extension
6057 Battery back-up facility
6054 Telescopic IV pole



Paper Roll Holder Fitted as Standard



Colour Olive

40 41



# **CLINNOVA** Gynae Pro - Models Available



#### **NV8483-CLS**

Electric Height Elevation & Tilt Electric Head Section Independently Braking Castors Classic Base



#### **NV8583-CLS**

Electric Height Elevation & Tilt Electric Head Section Retractable Castors Classic Base



#### **NV8683-CLS**

Electric Height Elevation & Tilt Electric Head Section Central Locking Castors Classic Base



#### **NV8483-PRM**

Electric Height Elevation & Tilt Electric Head Section Independently Braking Castors Premium Base



#### **NV8583-PRM**

Electric Height Elevation & Tilt Electric Head Section Retractable Castors Premium Base



#### **NV8683-PRM**

Electric Height Elevation & Tilt Electric Head Section Central Locking Castors Premium Base

42 43