



CLINNOVA Clinical 2 Section



A very strong and rigid couch available with a wide range of accessories. Suitable for supporting very heavy patients for examination and minor treatment procedures.

Models

See Page 6-7 for available models



Model NV2570-PRM
Colour Canard

Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted or electric backrest adjustable from -20° to +80°
- Anti-trap gas assisted backrest prevents risk of crushing and harm to patient, user or furniture
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working heights
- Linked retractable wheels with adjustable foot for uneven floors
- Lock & Go transport mode for easy moving
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Electric models include a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- 6004** Breathing hole and plug
- 9001** Perimeter foot switch (*NV2560 only*)
- 6200** Single foot switch (*NV2560 only*)
- 6207** Double foot switch (*NV2570 only*)
- 9002** Side support rails (Pair)
- 9002P** Side support rail cushions (Pair)
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6057** Battery back-up facility
- 6054** Telescopic IV pole



Paper Roll Holder
Fitted as Standard



Model NV2550-CLS
Colour Dark Blue

CLINNOVA Clinical 2 Section - Models Available



NV2550-CLS

Hydraulic Height Elevation
Gas Assisted Backrest
Classic Base



NV2560-CLS

Electric Height Elevation
Gas Assisted Backrest
Classic Base



NV2570-CLS

Electric Height Elevation
Electric Backrest
Classic Base



NV2550-PRM

Hydraulic Height Elevation
Gas Assisted Backrest
Premium Base



NV2560-PRM

Electric Height Elevation
Gas Assisted Backrest
Premium Base



NV2570-PRM

Electric Height Elevation
Electric Backrest
Premium Base