



CLINNOVA Clinical 3 Section



A very strong and rigid couch, with an additional adjustable section which increases versatility. Available with a wide range of accessories, suitable for supporting very heavy patients for examination and minor treatment procedures.

Models

See Page 10-11 for available models



Model NV3570-PRM
Colour Pepper

Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted or electric head section, adjustable from -20° to +80°
- Gas assisted or electric foot section, adjustable from 0° to +70°
- Anti-trap gas assisted sections prevent risk of crushing and harm to patient, user or furniture
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working heights
- Linked retractable wheels with adjustable foot for uneven floors
- Lock & Go transport mode for easy moving
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Electric models include a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colour

Optional Accessories

- 6004** Breathing hole and plug
- 9001** Perimeter foot switch (NV3560 only)
- 6200** Single foot switch (NV3560 only)
- 6209** Triple foot switch (NV3570 only)
- 9003** Side support rails (Pair)
- 9003P** Side support rail cushions (Pair)
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6057** Battery back-up facility
- 6054** Telescopic IV pole



Paper Roll Holder
Fitted as Standard



Model NV3550-CLS
Colour Burgundy

CLINNOVA Clinical 3 Section - Models Available



NV3550-CLS

Hydraulic Height Elevation
Gas Assisted Head Section
Gas Assisted Foot Section
Classic Base



NV3560-CLS

Electric Height Elevation
Gas Assisted Head Section
Gas Assisted Foot Section
Classic Base



NV3570-CLS

Electric Height Elevation
Electric Head Section
Electric Foot Section
Classic Base



NV3550-PRM

Hydraulic Height Elevation
Gas Assisted Head Section
Gas Assisted Foot Section
Premium Base



NV3560-PRM

Electric Height Elevation
Gas Assisted Head Section
Gas Assisted Foot Section
Premium Base



NV3570-PRM

Electric Height Elevation
Electric Head Section
Electric Foot Section
Premium Base