



CLINNOVA Mobile 2 Section



A very strong and rigid couch, fitted with large castors which make movement much easier. Fitted with a paper roll holder and side support rails as standard. Suitable for heavy patients for examination and minor treatment procedures.

Models

See Page 14-15 for available models



Model NV2670-PRM
Colour Dark Blue

Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted or electric backrest adjustable from -20° to +80°
- Anti-trap gas assisted backrest prevents risk of crushing and harm to patient, user or furniture
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working heights
- Choice of 125mm independent or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Electric models include a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- 6004** Breathing hole and plug
- 6200** Single foot switch (*NV2X60 only*)
- 6207** Double foot switch (*NV2X70 only*)
- 9002P** Side support rail cushions (Pair)
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6057** Battery back-up facility
- 6054** Telescopic IV pole



Paper Roll Holder
Fitted as Standard



Model NV2450-CLS
Colour Lotus Green

CLINNOVA Mobile 2 Section - Models Available



NV2450-CLS

Hydraulic Height Elevation
Gas Assisted Backrest
Independently Braking Castors
Classic Base



NV2460-CLS

Electric Height Elevation
Gas Assisted Backrest
Independently Braking Castors
Classic Base



NV2470-CLS

Electric Height Elevation
Electric Backrest
Independently Braking Castors
Classic Base



NV2650-CLS

Hydraulic Height Elevation
Gas Assisted Backrest
Central Locking Castors
Classic Base



NV2660-CLS

Electric Height Elevation
Gas Assisted Backrest
Central Locking Castors
Classic Base



NV2670-CLS

Electric Height Elevation
Electric Backrest
Central Locking Castors
Classic Base



NV2450-PRM

Hydraulic Height Elevation
Gas Assisted Backrest
Independently Braking Castors
Premium Base



NV2460-PRM

Electric Height Elevation
Gas Assisted Backrest
Independently Braking Castors
Premium Base



NV2470-PRM

Electric Height Elevation
Electric Backrest
Independently Braking Castors
Premium Base



NV2650-PRM

Hydraulic Height Elevation
Gas Assisted Backrest
Central Locking Castors
Premium Base



NV2660-PRM

Electric Height Elevation
Gas Assisted Backrest
Central Locking Castors
Premium Base



NV2670-PRM

Electric Height Elevation
Electric Backrest
Central Locking Castors
Premium Base